




PathSensors' Bio-Scan™ Assays use CANARY® Biosensor Technology



BioSensors Available

- | | |
|----------------------------------|--|
| <i>Salmonella spp.</i> | <i>Dengue virus</i> |
| <i>Listeria</i> | <i>RVF virus</i> |
| <i>E. coli</i> O157:H7 | <i>Vibrio cholerae</i> |
| <i>Bacillus anthracis</i> spores | (strains O139 & O1) |
| <i>Francisella tularensis</i> | <i>Chlamydia spp.</i> |
| <i>Yersinia pestis</i> | <i>Methicillin – Resistant Staphylococcus aureus</i> |
| <i>Orthopox virus</i> | Group B Strep |
| <i>Ricin toxin</i> | VEE virus |
| <i>Botulinum toxin</i> | Prion protein |
| <i>Brucella spp.</i> | FMD virus |
| <i>Bacillus subtilis</i> spores | <i>Shigella dysenteriae</i> |
| <i>Campylobacter</i> | - Digitonin labeled probes |
| <i>Ralstonia spp.</i> | - Double stranded DNA complex |
| <i>Potyvirus</i> | |
| <i>Phytophthora spp.</i> | |

NOTE: Not all biosensors are available on all PathSensors hardware platforms. Contact us for more information.

Bio-Scan™ Ricin Assay Kit

THE BIO-SCAN™ RICIN ASSAY KIT IS PART OF THE PATHSENSORS FAMILY OF BIO-SCAN™ ASSAY SOLUTIONS, WHICH OFFER HIGH SPEED, HIGH-SENSITIVITY PATHOGEN DETECTION.

The Bio-Scan™ Ricin Assay Kit allows for unparalleled detection of low levels of Ricin in less than 15 minutes. The Bio-Scan™ Ricin Assay Kit can be used with PathSensors' high performance Zephyr™ analysis platforms for a comprehensive and cost effective detection solution.

PathSensors' Bio-Scan™ Assays use breakthrough CANARY® technology originally developed by MIT-Lincoln Laboratory. Independent tests have shown that PathSensors' CANARY® provides superior speed to detection and sensitivity compared to technologies such as PCR and lateral flow. The technology utilizes biosensors created from B lymphocytes (white blood cells) modified to express surface-bound, target-specific antibodies and a bioluminescent protein. When the biosensor binds to its target pathogen, the antibodies trigger the intracellular release of calcium. This calcium causes the bioluminescent protein to emit light. Sophisticated algorithms analyze this light output, resulting in definitive "positive" or "negative" test results.

Advantages of this system are its extreme speed and sensitivity. The speed of detection is a result of rapid intracellular signaling. The sensitivity is achieved through signal amplification within the cell. The technology identifies targets in less than 12 minutes with analytical sensitivities down to picograms of target per sample.

Assay Specifications

Analytical Sensitivity (LoD)	400pg Ricin (ACAII, Ricinus communis agglutinin 60, RCA60, Castor Bean toxin)
Assay Sensitivity	1ng/mL Ricin
Resistance to Interferents	0 False Positives with 27 powders
Time to Results	12 minutes

NOTE: Performance validated by Battelle, Columbus, OH and Pacific Northwest National Laboratory, Richland, WA



Ricin Assay Protocol:

Dilute collected sample
in Sample Diluent



Add Ricin Capture Beads
and incubate with sample
(10 minutes)



Start Zephyr instrument,
enter operator ID &
assay barcode



Wash capture beads



Add Sample Diluent



Add Biosensors



Centrifuge 5 seconds



Transfer sample
to luminometer
(read for 2 minutes)



RETRIEVE RESULTS

Testing and Results:

The **Bio-Scan™ Ricin Assay** has been extensively tested for signal specificity with a variety of suspicious powders.

The suspicious powders used for testing were:

Acetaminophen	Borax®	Epsom salt	MiraLAX®	Tooth powder
Aerosil powder	Chalk	Foot powder	Powdered sugar	Tums™
Ajax®	Corn starch	Gym chalk	Road dust	Yeast
Baking powder	Dairy creamer	Infant formula	Salt powder	
Baking soda	DiPel	Instant pectin	Spackling powder	
Bentotite	Dry milk	Kaolin	Talcum powder	

SOURCE: Suspicious Powders Panel, Critical Reagents Program.

RICIN ASSAYS		
Supplier: Device	Approx. Time (min)	Reported LoD (ng/mL)
CANARY		
PathSensors: Zephyr™/Bio-Scan™	< 12	1 ng
IMMUNOASSAY		
Typical Immunoassay Device	5 - 20	Typical 1-10 ng
PCR		
Typical PCR Device	30-60	Cannot directly detect

* Table generated using information found in the Biodetection Technologies for First Responders: 2014 Edition; RM Ozanich, CL Baird, RA Barholomew, HA Colburn, TM Straub, CJ Bruckner-Lea, Pacific Northwest National Laboratory. <http://biodetectionresource.pnnl.gov>

For copies of independent tests, contact PathSensors, Inc.

About PathSensors:

PathSensors is a leading biotechnology solutions and environmental testing company, providing high speed, high sensitivity, pathogen and threat detection solutions.

PathSensors' solutions can detect a wide range of threats, including anthrax, ricin, *Ebola* and *salmonella*. PathSensors' technology is being used today by government and commercial customers for multiple applications.

